

Rodin® Denture Base

Denture Base Resin

Rodin® Denture Base Resin is designed to be used alongside Rodin® Sculpture and is suitable for printing all types of denture bases. Maximizes durability to match performance of traditional denture acrylics, outperforming all other non-filled denture base resins that may otherwise break upon dropping.

Benefits & Features



Unmatched Strength



Highest Fracture & Wear Resistance



Matches Color of Gingiva



Stunning Aesthetics



FDA 510K Cleared & Biocompatible



High translucency on Denture Base and Sculpture. Courtesy of Lee Culp, CDT, Sculpture Studios



Pink shade mimics gingiva tissue. Courtesy of Daniel Nowak, CDT

Technical Data

	Rodin® Denture Base
Viscosity 23C (cps)	800 - 1200 cP
Shore Hardness ASTM D2240	≥ 95 Shore D
Flex Strength MPa (SD) ISO 20795-1	93 - 98 MPa
Flexural Modulus (MPa)	2038 - 2573 MPa

Order Information

3DR-DB-ORF

Rodin® Denture Base Original, 1 kg Bottle

Rodin® Resin

ADDITIONAL PRODUCT LINE



3DR-TRAYF

Rodin® Bite Tray

Custom Bite Tray Resin, 1 kg



3DR-SGF

Rodin® Surgical Guide

Surgical Guide Resin, 1 kg



3DR-TANF

Rodin® Model

Model Resin, 1 kg



3DR-TIF

Rodin® Try-In

Try-in Denture Resin, 1 kg



3DR-OIBTF

Rodin® Ortho IBT

Indirect Bonding Tray Resin, 1 kg



3DR-CASTF

Rodin® Castable

Castable Resin, 1 kg

V

Authorized Reseller



Aurident
Precious smiles for life.

aurident.com
800-422.7373