

Rodin[®] Splint

Clear Bite Splint Resin

Rodin[®] Splint is designed for printing a variety of rigid dental splints, night guards, and retainers. Due to its very high strength and resistance, it is ideal for patients suffering from acute bruxism and night grinding and can be prescribed for both bite splints and snoring appliances.

Benefits & Features



Unmatched Strength



Highest Fracture & Wear Resistance



Minimal Yellow Discoloration



Wide Range of Applications



FDA 510K Cleared & Biocompatible



Finished print of Rodin Splint
Courtesy of Rick Ferguson, DMD

Technical Data

	Rodin[®] Splint
Viscosity 23C (cps)	400 - 650 cP
Shore Hardness ASTM D2240	≥ 95 Shore D
Flex Strength MPa (SD) ISO 20795-1	80 - 90 MPa
Flexural Modulus (MPa)	3527 - 4228 MPa

Order Information

3DR-SPLINTF Rodin[®] Splint, 1 kg Bottle

Rodin[®] Resin

ADDITIONAL PRODUCT LINE



3DR-TRAYF

Rodin[®] Bite Tray

Custom Bite Tray Resin, 1 kg



3DR-SGF

Rodin[®] Surgical Guide

Surgical Guide Resin, 1 kg



3DR-TANF

Rodin[®] Model

Model Resin, 1 kg



3DR-TIF

Rodin[®] Try-In

Try-in Denture Resin, 1 kg



3DR-OIBTF

Rodin[®] Ortho IBT

Indirect Bonding Tray Resin, 1 kg



3DR-CASTF

Rodin[®] Castable

Castable Resin, 1 kg

V

Authorized Reseller



Aurident

Precious smiles for life.

aurident.com

800-422.7373