



Rodin® SOFT RELINE

Compatible with 3D, PMMA & Acrylic Denture

NEW

Rodin® Soft Reline features a unique silicone-based formulation that improves the fit and enhances the performance of both removable and fixed full or partial dentures. Tailored for superior bonding and tear strength, it is the ideal reline solution for 3D printed polymer bases, as well as PMMA milled and traditional acrylic materials.

Benefits & Features

- Exceptional Adhesion and Tear Strength
- Ease of Workflow
- Superior Performance
- Esthetic and Patient Comfort

REF 24006 Soft Reliner Kit

Components




VPS Adhesive Primer & VPS Release

VPS Adhesive boasts superior adhesion, thanks to its innovative formulation that penetrates deeply into denture materials to form a resilient bond. Meanwhile, VPS Adhesive Release dissolves tough adhesives without residue, ensuring a comfortable and hassle-free cleaning experience for users.

10mL VPS Adhesive Primer **REF** 24008

10mL VPS Release **REF** 24009



VPS Reline Cartridge

Unique silicone-based reline formulation of VPS Reline Cartridge improves the fit and enhances the performance of both removable and fixed full or partial dentures. Its physical properties make it the ideal solution for 3D-printed polymer bases, as well as PMMA milled and traditional acrylic materials.

REF 24007 50mL Reline Silicone Cartridge

REF 240013 Wide Mouth Mixing Tips Refil, 48/pk



Polishing Brush & Tapered Stone

Specialized polishing tools are included to complement the needs of denture maintenance. The polishing brush assists in refining the reline surface to a flawless finish, while the tapered stone bur shapes materials with precision.

REF 240010 Tapered Cone Finishing Stone

REF 240011 Brown Polishing Wheel