

 RODIN®  
**CHROMA FLASH™**

*Consistent Flash Post-Curing Designed with Rodin® Materials*



Stable Material  
Performance



Predictable Surface  
and Aesthetic Outcomes



Open, Material-First  
System



Simplified, Low-Sensitivity  
Workflow

Engineered to deliver consistent, reliable polymerization without nitrogen, Chroma Flash simplifies post-curing while supporting predictable mechanical properties, surface quality, and workflow outcomes.



## Comparison Table

|                                | Chroma Flash | Leading competitors |
|--------------------------------|--------------|---------------------|
| Xenon Flash Curing             | ✓            | X                   |
| Integrated Material Library    | ✓            | X                   |
| Customizable Save Profiles     | ✓            | X                   |
| Tack Cure Cycles               | ✓            | X                   |
| Seamless Workflow Intergration | ✓            | X                   |
| Wireless Connectivity          | ✓            | X                   |
| N2 Compatibility               | ✓            | X                   |

## Technical Specifications

| CATEGORY                       | SPECIFICATIONS        |
|--------------------------------|-----------------------|
| Inner Curing Chamber Dimension | 111 x 111 x 43.5 mm   |
| Device Dimensions              | 256 x 330 x 223.76 mm |
| Light / Source                 | Flash / Stroboscopic  |
| Wavelength                     | 280-950 nm            |
| Rated Voltage                  | 110-240V              |
| Weight                         | 8.2 kg                |
| Open / Close System            | Open                  |
| N <sub>2</sub> Compatibility   | Yes                   |



DESIGN



PRINT



MATERIALS



WASH & FINISH



CURE



DELIVER

